



The **bionivid** Science Blog

WGS LIBRARY PREPARATION: A BRIEF OVERVIEW

JULY - 2025 - VII



- Team Genomics Lab, Bionivid


BIONIVID TECHNOLOGY PRIVATE LIMITED
BENGALURU, KARNATAKA, INDIA.

✉ info@bionivid.com | sales@bionivid.com

☎ +91 95356 19191 | +91 96064 56771



Whole Genome Sequencing (WGS) is a transformative technology in genomics, providing an unbiased, high-resolution view of the complete DNA sequence of an organism. Among the various formats of WGS, **short-insert library preparation paired with 150×2 paired-end sequencing chemistry** has become a gold standard in both research and clinical genomics. This method strikes an effective balance between throughput, accuracy, and cost-efficiency, enabling broad adoption across a wide range of biological investigations.



What is Short-Insert WGS?

The Short-insert WGS refers to the sequencing of DNA fragments that are typically around 250–300 base pairs (bp) in length, not including adapter sequences. When these libraries are sequenced using paired-end 150 bp reads, both ends of each DNA fragment are sequenced, resulting in 150×2 read pairs. This paired-end strategy significantly improves the ability to accurately map reads to a reference genome, resolve repetitive regions, and detect small-scale genomic variations such as single nucleotide polymorphisms (SNPs) and short insertions or deletions (indels).



Enzymatic Fragmentation: A Controlled Approach

Traditionally, mechanical methods like ultrasonication (e.g., Covaris) have been used to fragment DNA. However, enzymatic fragmentation is increasingly favored for its convenience, reproducibility, and compatibility with automation.

Enzymatic methods utilize specialized DNA fragmentase enzymes that cleave genomic DNA at random locations under controlled reaction conditions. This approach eliminates the need for expensive instrumentation, reduces sample handling, and is gentler on the DNA, an advantage when working with limited or sensitive samples.

The degree of fragmentation can be fine-tuned by adjusting enzyme concentration, reaction time, and temperature, allowing the generation of libraries with a precise insert size distribution.

Importantly, enzymatic fragmentation can be integrated into streamlined library preparation workflows, including single-tube protocols, which minimize DNA loss and contamination risk.



Library Preparation Workflow

Once DNA is fragmented enzymatically, it undergoes end-repair and A-tailing to prepare the blunt or overhanging ends for ligation. This step ensures that all fragments have a compatible 3' adenine overhang, which is essential for ligation to T-tailed sequencing adapters.


Next, adapter ligation is performed, attaching platform-specific adapters (e.g., for Illumina systems) to both ends of each DNA fragment. These adapters contain necessary elements for hybridization to the flow cell, sequencing priming sites, and often sample-specific barcodes (indexes) that allow for multiplexing. The ligation reaction is followed by a bead-based purification step (using AMPure XP or equivalent) to remove free adapters and very short DNA fragments.

To generate sufficient material for sequencing, the adapter-ligated library is subjected to limited-cycle PCR enrichment. The number of cycles is minimized (typically 5–7) to avoid introducing bias or over-amplification artifacts.

The PCR enrichment reaction is followed by a bead-based purification step (using AMPure XP or equivalent) to remove any PCR artifacts. The PCR enriched library can now be taken up for Library QC.



Library Preparation Workflow



Quality Control and Library Assessment

A critical aspect of short-insert WGS library preparation is thorough quality control (QC). Libraries are quantified using fluorometric assays such as Qubit dsDNA HS, which specifically measures double-stranded DNA. Size distribution and integrity are evaluated using capillary electrophoresis platforms like Agilent TapeStation or Bioanalyzer. A successful short-insert WGS library will typically show a tight size distribution centered around 400–450 bp, accounting for both insert and adapter sequences.

Libraries must meet both concentration and size criteria to proceed to sequencing. Undersized or low-concentration libraries may lead to inefficient cluster formation or low data yield. In some workflows, additional normalization or bead-based size selection may be performed to optimize performance.

Library preparation is a critical foundation for the success of any NGS-based experiment, particularly WGS, directly influencing data quality, coverage uniformity, and variant detection accuracy. As the bridge between raw genomic DNA and high-throughput sequencing, this step determines how truly the genome is represented in the final data. Precision in fragmentation, adapter ligation, and amplification ensures that the complexity and integrity of the genome are preserved. Without a well-optimized and quality-controlled library, even the most advanced sequencing platforms cannot deliver meaningful or interpretable results.



bionivid

BIONIVID TECHNOLOGY PRIVATE LIMITED

BENGALURU, KARNATAKA, INDIA.

✉ info@bionivid.com | sales@bionivid.com

☎ +91 95356 19191 | +91 96064 56771



www.bionivid.in



The Genome Education
Nurturing Careers